

Dyna Cone Plus⁺

Dyna Cone Plus⁺ Elite

This design is for cones that are centrally positioned and have a smaller diameter with varying optical zones. Start with an 8.8mm diameter, but 8.0, 8.4, and 9.2 are also available. The peripheral system is a heavily blended series of spherical cuts for better centration and maximum alignment. The Dyna Cone Plus+ and Dyna Cone Plus+ Elite have a 0.4mm larger optical zone than the Dyna cone.

Features	Benefits
Cone Contour Technology	More closely matches the back surface to the cone profile
Spherical cuts, heavily blended series of cuts	Controlled positioning and reproducibility
Can be made steeper, up to 4.50mm	Flexibility; ability to fit any centrally located nipple cone
Edge lifts - standard, flat, steep	Better alignment in the periphery
Proven peripheral design system	Fewer edge lift adjustments
Front toric and bitoric options	Fits the most complex and challenging cases (corrects residual astigmatism)
Quad Sym availability	Resolves 6 o'clock standoff and pop-off
Diameters available 8.0 to 9.9	Fits a wide range of patient fissures
Available in all GP materials	Fitter material preference (Elite features Boston XO)
Professional Consultation team	Reduced chair time and increased success
Quick turnaround time	Patients treated/fit faster
Loaner trial sets available	Efficient, first time fitting success



Kansas City, Missouri
 (800) 228-2691 • (303) 237-6927 local
 Fax: (800) 661-6707 • (303) 274-6707 local
orderentry@lensdynamics.com

