

Dyna Semi Scleral

This lens comes in diameters 13.5, 14.0, 14.5, and 15.0mm. Standard BOZ is 9.0mm but may be changed. Edge lifts may be increased or decreased as needed. +10.00 to -20.00 powers

Features	Benefits
Scleral fit/no corneal touch	Provides soft lens comfort
Larger diameters available 13.5 to 15.0	Better comfort and positioning
Edge lifts - standard, flat, steep	Better alignment in the periphery
Proven peripheral design system	Fewer edge lift adjustments
Available in a front aspheric multifocal	Presbyopic correction
Reverse geometry option	Allows for a closer fit for irregular corneal profile
Single lens material	Avoids hybrid lens difficulties
Available in all high dk GP materials	Fitter preference of high dk materials
Professional Consultation team	Reduced chair time and increased success
Quick turnaround time	Patients treated/fit faster
Loaner trial sets available	Efficient, first time fitting success



Innovators in Irregular Corneal Treatment

Kansas City, Missouri
(800) 228-2691 • (303) 237-6927 local
Fax: (800) 661-6707 • (303) 274-6707 local
orderentry@lensdynamics.com

