

LENS DYNAMICS

PRECISION LENSES

Innovators in Irregular Corneal Treatment

Kansas City, Missouri

Phone: (800) 228-2691 ♦ (303) 237-6927 **local**

Fax: (800) 661-6707 ♦ (303) 274-6707 local
reception@lensdynamics.com

FIRESTONE OPTICS, INC.

The Standard of Excellence



Dyna Irregular Cornea Lens Standard Trial Set Programs

Loaner Set	<ul style="list-style-type: none"> We send the set to you; you are billed for the set when it goes out. You fit your patient Return the set to us and you are credited back for the billing less any shipping charges.
Earned Set	<ul style="list-style-type: none"> We send the set to you; you are billed for the set but given 90 day terms. Fit 3 successful patients (6 lenses) within the 90 day billing period. Request a credit within the 90 days and you will never be billed for the set.
Purchased Set	We have many set sizes available at reasonable prices to meet your every need.

Design	Diameter <i>Choose One</i>	Series	# of Lenses	Base Curves	Price
Dyna Cone® Plus Dyna Cone Plus Elite	8.0, 8.4, 8.8, 9.2, 9.6, 9.9	Standard	26	5.10, 5.20, 5.30, 5.40, 5.50, 5.60, 5.70, 5.80, 5.90, 6.00, 6.10, 6.20, 6.30, 6.40, 6.50, 6.60, 6.70, 6.80, 6.90, 7.00, 7.10, 7.20, 7.30, 7.40, 7.50, 7.60	\$299.00
		Flatter	5	7.60, 7.70, 7.80, 7.90, 8.00	\$57.50
	8.0, 8.4, 8.6	Steeper	6	4.50, 4.60, 4.70, 4.80, 4.90, 5.00	\$69.00
	8.8	Steeper	5	4.60, 4.70, 4.80, 4.90, 5.00	\$57.50
	9.2, 9.4, 9.6, 9.8	Steeper	4	4.70, 4.80, 4.90, 5.00	\$46.00
Dyna Intra-Limbal® Dyna Intra-Limbal Elite	10.8, 11.2, 11.6, 12.0	Standard	14	6.49, 6.62, 6.75, 6.89, 7.03, 7.18, 7.34, 7.50, 7.67, 7.85, 8.04, 8.23, 8.44, 8.65	\$161.00
	10.8, 11.2, 11.6	Flatter	6	8.77, 8.88, 9.00, 9.12, 9.25, 9.38	\$69.00
	10.8	Steeper	6	5.82, 5.92, 6.03, 6.14, 6.25, 6.37	\$69.00
	11.2	Steeper	5	5.92, 6.03, 6.14, 6.25, 6.37	\$57.50
	12.0	Steeper	4	6.03, 6.14, 6.25, 6.37	\$46.00
Dyna Intra-Limbal® <i>special reverse geometry</i>	10.8, 11.2, 11.6	Standard	18 Series A, B, C 18 Series E, F, G	8.04, 8.23, 8.44, 8.65, 8.88, 9.12	\$207.00
	10.8, 11.2, 11.6	Flatter	18 Series A, B, C 18 Series E, F, G	9.38, 9.64, 9.93, 10.23, 10.55, 10.89	\$207.00
	10.8, 11.2, 11.6	Steeper	18 Series A, B, C 18 Series E, F, G	7.03, 7.18, 7.34, 7.50, 7.67, 7.85	\$207.00
Dyna Semi-Scleral™	13.5, 14.0, 14.5, 15.0, 15.8	Standard	13	5.50, 5.70, 5.90, 6.10, 6.30, 6.50, 6.70, 6.90, 7.10, 7.30, 7.50, 7.70, 7.90	\$199.00
Post Graft Floating OZ	10.4 - floating OZ	Standard	12	6.92, 7.06, 7.21, 7.37, 7.53, 7.70, 7.88, 8.07, 8.26, 8.47, 8.68, 8.91	\$138.00
	10.4 - 9.0 OZ	Standard	16	6.37, 6.49, 6.61, 6.75, 6.89, 7.03, 7.18, 7.34, 7.50, 7.67, 7.85, 8.04, 8.23, 8.44, 8.65, 8.88	\$184.00
Post Lasik	10.0 with 8.0 OZ	Standard	18 Series A, B, C 18 Series E, F, G	8.04, 8.23, 8.44, 8.65, 8.88, 9.12	\$207.00